



AccuRead Arm DB61M (Blood Pressure Monitor & WHO Light Sign Indicator)

TRULY INSTRUMENT

Product Qualification

CE0434: This product fulfils the requirement for Directive 93/42/EEC concerning Medical Device.

FDA KO14141/ISO 9001/MDD Annex V:

Manufactured in accordance with International standards

Specifications

Blood Pressure

Measurement Range: 20~280mm Hg
Accuracy: +/- 3mmHg
Display: Systole, Diastole, Pulse


Pulse

Measurement Range: 40~195 pulse/min.
Pulse Rate Accuracy: +/- 5% of the reading

Power Dimension

Power: 3V, 4 x AAA size
Dimension: 4.8" (W) x 6" (H) x 5.3" (D) Weight: Approx. 344g (excl. batteries), cuff 120g
Arm Circumference : 22 - 34 cm (Special order for extra size)
Auto Power off after 3 minutes

PC-Link: This unit also offered PC-link software which include USB cable and software cd.



Light Sign Indicator Chart